

PRODUCT INFORMATION

Target	DLL3
Synonyms	SCDO1
Description	Recombinant Cynomolgus DLL3(426-488) protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	A0A2K5WSR1
Expression Host	HEK293
Tag	C-Human Fc tag
Molecular Characterization	DLL3(Glu426-Arg488)+hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 32.6 kDa after removal of the signal peptide. The apparent molecular mass of cDLL3(426-488)-hFc is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Cynomolgus DLL3(426-488) Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

