

PRODUCT INFORMATION

Tag	C-Flag Tag
Target	AGRF5
Synonyms	GPR116, KPG_001
Description	Human AGRF5 full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q8IZF2
Expression Host	HEK293
Protein Families	GPCR,Transmembrane,Druggable Genome,
Protein Pathways	GPCRDB Other,
Molecular Weight	The human full length AGRF5 protein has a MW of 149.5kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Receptor that plays a critical role in lung surfactant homeostasis. May play a role in controlling adipocyte function.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

