

PRODUCT INFORMATION

Target	ANGPTL7
Synonyms	AngX;CDT6:dj647M16.1
Description	Recombinant Human ANGPTL7 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	O43827
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	ANGPTL7(Gln27-Pro346) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 37.9 kDa after removal of the signal peptide. The apparent molecular mass of ANGPTL7-His is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	Has a role in the formation and organization of the extracellular matrix. In the eye, it functions as a mediator of dexamethasone-induced matrix deposition in the trabecular meshwork, the tissue responsible for the outflow of the ocular aqueous humor and for the maintenance of intraocular pressure (PubMed:21199193). Is a negative regulator of angiogenesis in the cornea, and plays a major role in maintaining corneal avascularity and transparency (PubMed:25622036).[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated





Figure 1. Human ANGPTL7 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

