

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	CNGA2
Synonyms	CNCA, CNCA1, CNG2, OCNC1, OCNCALPHA, OCNCa
Description	Human CNGA2-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q16280
Expression Host	HEK293
Protein Families	Ion Channels: Cyclic nucleotide gated
Protein Pathways	N/A
Molecular Weight	The human full length CNGA2-Strep protein has a MW of 76 kDa Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Formulation & Reconstitution	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Storage&Shipping	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Sterility	The protein encoded by this gene represents the alpha subunit of a cyclic nucleotide-gated olfactory channel. The encoded protein contains a carboxy-terminal leucine zipper that mediates channel formation. [provided by RefSeq, Jan 2010]
Background	
Usage	Research use only
Conjugate	Unconjugated

