

PRODUCT INFORMATION

Tag	C-Flag Tag
Target	CTSR4
Synonyms	N/A
Description	Human CTSR4 full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q7RTX7
Expression Host	HEK293
Protein Families	Ion Channels: Other
Protein Pathways	N/A
Molecular Weight	The human full length CTSR4 protein has a MW of 54.1kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Voltage-gated calcium channel that plays a central role in calcium-dependent physiological responses essential for successful fertilization, such as sperm hyperactivation, acrosome reaction and chemotaxis towards the oocyte.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

