

## PRODUCT INFORMATION

<b>Target</b>	IL-28B
<b>Synonyms</b>	Interleukin-28B;IL-28B;Cytokine Zcyto22;Interferon Lambda-3;IFN-Lambda-3;Interferon Lambda-4;IFN-Lambda-4;Interleukin-28C;IL-28C;IL28B;IFNL3;IFNL4;IL28C;ZCYTO22
<b>Description</b>	Recombinant Human Interleukin-28B is produced by our Mammalian expression system and the target gene encoding Val22-Val196 is expressed.
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	Q8IZI9
<b>Expression Host</b>	HEK293
<b>Tag</b>	
<b>Molecular Characterization</b>	Not available
<b>Molecular Weight</b>	19.6 KDa
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Sterility</b>	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
<b>Background</b>	Interleukin-28B, also known as Cytokine Zcyto22, Interferon lambda-3, Interferon lambda-4, IFNL3, IFNL4, ZCYTO22 and IL28B, is a secreted cytokine which belongs to the IL-28/IL-29 family. IL-28 has also been shown to play a role in the adaptive immune response. IL28B has immunomodulatory activity and up-regulates MHC class I antigen expression. IL28B displays potent antiviral activity and antitumor activity. In addition, IL28B is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

DIMABIO CONFIDENTIAL

