

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	MCLN3
Synonyms	TRP-ML3, TRPML3
Description	Human MCLN3-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q8TDD5
Expression Host	HEK293
Protein Families	Ion Channels: Transient receptor potential
Protein Pathways	N/A
Molecular Weight	The human full length MCLN3-Strep protein has a MW of 64.2 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	This gene encodes one of members of the mucolipin cation channel proteins. Mutation studies of the highly similar protein in mice have shown that the protein is found in cochlea hair cells, and mutant mice show early-onset hearing loss and balance problems. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]
Usage	Research use only
Conjugate	Unconjugated

