

PRODUCT INFORMATION

Target	SCTR
Synonyms	SR; SECR
Description	Recombinant human SCTR Protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P47872
Expression Host	HEK293
Tag	C-Human Fc tag
Molecular Characterization	SCTR(Leu28-Lys141) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 39.2 kDa after removal of the signal peptide.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The protein encoded by this gene is a G protein-coupled receptor and belongs to the glucagon-VIP-secretin receptor family. It binds secretin which is the most potent regulator of pancreatic bicarbonate, electrolyte and volume secretion. Secretin and its receptor are suggested to be involved in pancreatic cancer and autism. [provided by RefSeq, Jul 2008]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human SCTR Protein, hFc Tag on SDS-PAGE under reducing condition.

