

**PRODUCT INFORMATION**

<b>Tag</b>	C-Flag Tag
<b>Target</b>	SMO
<b>Synonyms</b>	CRJS, FZD11, Gx, PHLS, SMOH
<b>Description</b>	Human SMO full length protein-synthetic nanodisc
<b>Delivery</b>	6~8weeks
<b>Uniprot ID</b>	Q99835
<b>Expression Host</b>	HEK293
<b>Protein Families</b>	ES Cell Differentiation/IPS ,Transmembrane,Stem cell relevant signaling - DSL/Notch pathway,Druggable Genome,
<b>Protein Pathways</b>	GPCRDB Other,Hedgehog Netpath 10,Hedgehog Netpath 10,Apoptosis,Cancer,Notch,
<b>Molecular Weight</b>	The human full length SMO protein has a MW of 86.4kDa
<b>Formulation &amp; Reconstitution</b>	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	The protein encoded by this gene is a G protein-coupled receptor that interacts with the patched protein, a receptor for hedgehog proteins. The encoded protein transduces signals to other proteins after activation by a hedgehog protein/patched protein complex. [provided by RefSeq, Jul 2010]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

