

PRODUCT INFORMATION

| | |
|---|--|
| Tag | C-Flag Tag |
| Expression Host | HEK293 |
| Target | TAAR2 |
| Synonyms | GPR58, taR-2 |
| Description | Human TAAR2 full length protein-synthetic nanodisc |
| Uniprot ID | Q9P1P5 |
| Protein Families | GPCR,Transmembrane,Druggable Genome, |
| Protein Pathways | GPCRDB Other, |
| Molecular Weight | The human full length TAAR2 protein has a MW of 40.1kDa |
| Delivery | 6~8weeks |
| Formulation & Reconstitution | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution |
| Sterility | Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use. Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Storage&Shipping | |
| Background | Orphan receptor.[UniProtKB/Swiss-Prot Function] |
| Usage | Research use only |
| Conjugate | Unconjugated |

