

PRODUCT INFORMATION

Target	TFRC
Synonyms	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46
Description	Recombinant human TFRC(121-760) Protein with N-terminal 10×His tag
Delivery	In Stock
Uniprot ID	P02786
Expression Host	HEK293
Tag	N-10×His tag
Molecular Characterization	10×His tag TFRC(Arg121-Phe760)
Molecular Weight	The protein has a predicted molecular mass of 73.1 kDa after removal of the signal peptide.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	This gene encodes a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Multiple alternatively spliced variants have been identified. [provided by RefSeq, Sep 2015]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human TFRC(121-760) Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

