

PRODUCT INFORMATION

Target	CD45
Description	The embedded and sliced tumor sections (4µm section) from M-NSG mice implanted with Daudi cells
Host Cells	Daudi
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage&Shipping	Room temperature
Usage	Research use only




Figure 1. IHC analysis of M-NSG Daudi DiSliceX™ SlideSet section (Cat. No. SLI100008) with Anti-CD38 antibody(DM28)Rabbit mAb (Cat. No. DME100027) at a concentration of 10 µg/ml. The inset image shows the absence of staining in secondary antibody-only control.



DIMABIO CONFIDENTIAL

