

PRODUCT INFORMATION

Target	B7H5, BST1, CCR2, CD33, CD70, CD93, CHI3L1, TREM2
Description	The embedded and sliced tumor sections (4μm section) from M-NSG mice implanted with THP-1 cells
Host Cells	THP-1
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only

