

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | PF-05082566, PF-2566 |
| Conjugate | PE-conjugated |
| Synonyms | TNFRSF9;4-1BB;CD137;CDw137;ILA |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | Human IgG2 - Lambda1 |
| Reactivity | Human |
| Target | 4-1BB |
| Uniprot ID | Q07011 |
| Description | PE-conjugated Anti-4-1BB (utomilumab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

