

PRODUCT INFORMATION

Common Name	H8, mAb5T4
Conjugate	PE-conjugated
Synonyms	M6P1;WAIF1;TPBG
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1 μ l/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	5T4
Uniprot ID	Q13641
Description	PE-conjugated Anti-5T4 (H8 biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-5T4 (H8 biosimilar) mAb on HEK293 cells (Blue histogram) or Jurkat cells (Red histogram).

DIMABIO CONFIDENTIAL

