

## PRODUCT INFORMATION

<b>Common Name</b>	IMGC 936,IMGC-936,IMGC936
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	CORD9;MCMP;MDC9;Mltng
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG1(K214R,M252Y,S254T,T256E,S442C) - Kappa
<b>Reactivity</b>	Human
<b>Target</b>	ADAM9
<b>Uniprot ID</b>	Q13443
<b>Description</b>	PE-conjugated Anti-ADAM9(izeltabart biosimilar) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

