

PRODUCT INFORMATION

Clone ID DM153 B7-H6 **Target**

B7-H6;NCR3LG1;B7 Homolog 6 **Synonyms**

Host Species Rabbit

PE-conjugated Anti-B7-H6 antibody(DM153); Description

Rabbit mAb

Delivery Under Development

Uniprot ID Q68D85 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Reconstitution

Background

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

B7-H6 belongs to the B7 family (see MIM 605402)

and is selectively expressed on tumor cells. Interaction of B7H6 with NKp30 (NCR3; MIM

611550) results in natural killer (NK) cell activation and cytotoxicity.

Usage Research use only

Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or

> > Email: info@dimabio.com Website: www.dimabio.com

DIMA Disclaimer reverse engineering attempt is prohibited. We are

actively scrutinizing all patent application to ensure no IP infringement.

