

**PRODUCT INFORMATION**

<b>Common Name</b>	LymphoStat-B, hBlySmAb-1.1
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	TNFSF13B;BAFF;BLYS;CD257;DTL;TALL1;THANK;TNFSF20;ZTNF4;TALL-1
<b>Applications</b>	Flow Cyt
<b>Endotoxin</b>	Less than 1.0 EU/ $\mu$ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Homo sapiens
<b>IgG type</b>	Human IgG1 - Lambda2
<b>Reactivity</b>	Human
<b>Target</b>	BAFF
<b>Uniprot ID</b>	Q9Y275
<b>Description</b>	PE-conjugated Anti-BAFF (belimumab biosimilar) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Sterility</b>	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 $\mu$ m) prior to use.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

