

## PRODUCT INFORMATION

<b>Common Name</b>	2857916, GSK2857916, J6M0, J6M0-mcMMAF
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	CD269;TNFRSF17;BCM;BCMA
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG1 - kappa
<b>Reactivity</b>	Human
<b>Target</b>	BCMA
<b>Uniprot ID</b>	Q02223
<b>Description</b>	PE-conjugated Anti-BCMA (belantamab biosimilar) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

