

PRODUCT INFORMATION

Common Name	WX-G250, cG250
Conjugate	PE-conjugated
Synonyms	CAIX;MN
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	CA9
Uniprot ID	Q16790
Description	PE-conjugated Anti-CA9 (girentuximab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

