

PRODUCT INFORMATION

Common Name	CAN34D3-65,PF-06443537
Conjugate	PE-conjugated
Synonyms	IL-31
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	IgG1
Reactivity	Canis lupus familiaris
Target	IL31
Uniprot ID	C7G0W1
Description	PE-conjugated Anti-Canine IL31(lokivetmab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

