Cat. No. DMC100477P



## **PRODUCT INFORMATION**

Clone ID **DMC477 Target** CCR6

BN-1; C-C CKR-6; CC-CKR-6; CCR-6; CD196; CKR-**Synonyms** 

L3; CKRL3; CMKBR6; DCR2; DRY6; GPR29;

GPRCY4; STRL22

**Host Species** Rabbit

PE-conjugated Anti-CCR6 antibody(DMC477); **Description** 

IgG1 Chimeric mAb Under Development

P51684 **Uniprot ID** 

Rabbit/Human Fc chimeric IgG1 IgG type

Monoclonal Clonality Reactivity Human **Applications** Flow Cyt

Recommended

**Background** 

**DIMA Disclaimer** 

Delivery

Flow Cyt 1:100 **Dilutions** 

Purified from cell culture supernatant by affinity **Purification** 

chromatography

Formulation & Liquid PBS with 0.05% Proclin300, 1% BSA Reconstitution

Store at 2°C-8°C for 6 months Storage & Shipping

> This gene encodes a member of the beta chemokine receptor family; which is predicted to be a seven transmembrane protein similar to G protein-coupled receptors. The gene is

preferentially expressed by immature dendritic cells and memory T cells. The ligand of this receptor is macrophage inflammatory protein 3 alpha (MIP-3 alpha). This receptor has been shown to be important for B-lineage maturation and antigen-driven B-cell differentiation; and it

may regulate the migration and recruitment of dentritic and T cells during inflammatory and immunological responses. Alternatively spliced transcript variants that encode the same protein have been described for this gene. [provided by

RefSeq; Jul 2008]

Usage Research use only

Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing apparent application to

ensure no IP infringement.





