

PRODUCT INFORMATION

Common Name	CSL-362-AML, CSL362, JNJ-56022473
Conjugate	PE-conjugated
Synonyms	CD123;IL3R;IL3RA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 – kappa
Reactivity	Human
Target	CD123
Uniprot ID	P26951
Description	PE-conjugated Anti-CD123 (talacotuzumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

