

PRODUCT INFORMATION

Common Name	IC14
Conjugate	PE-conjugated
Synonyms	CD14
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	Human IgG4(S228P) - Kappa
Reactivity	Human
Target	CD14
Uniprot ID	P08571
Description	PE-conjugated Anti-CD14(atibuclimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

