

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | RB4v1.2 |
| Conjugate | PE-conjugated |
| Synonyms | CD19,B4,CVID3,MGC12802 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Chimeric |
| IgG type | Human IgG1 - kappa |
| Reactivity | Human |
| Target | CD19 |
| Uniprot ID | P15391 |
| Description | PE-conjugated Anti-CD19(loncastuximab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

