

PRODUCT INFORMATION

Common Name	ICP-B04;ICP B04;MOR 208;MOR208;MOR-208;XENP-5574;XENP5574;XENP 5574;XmAb-5574;XmAb5574;XmAb 5574
Conjugate	PE-conjugated
Synonyms	CD19;B4;CVID3;MGC12802
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1 μ l/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	N/A
IgG type	Human IgG1-G2 - kappa
Reactivity	Human
Target	CD19
Uniprot ID	P15391
Description	PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb on Raji cells (Blue histogram) or 293T cells (Red histogram).

