

**PRODUCT INFORMATION**

<b>Common Name</b>	ICP-B04;ICP B04;MOR 208;MOR208;MOR-208;XENP-5574;XENP5574;XENP 5574;XmAb-5574;XmAb5574;XmAb 5574
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	CD19;B4;CVID3;MGC12802
<b>Applications</b>	Flow Cyt
<b>Endotoxin</b>	Less than 1.0 EU/ $\mu$ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	Flow Cyt 1 $\mu$ l/test
<b>Formulation &amp; Reconstitution</b>	Liquid,PBS with,0.18%BSA□0.1%Proclin300
<b>Host Species</b>	N/A
<b>IgG type</b>	Human IgG1-G2 - kappa
<b>Reactivity</b>	Human
<b>Target</b>	CD19
<b>Uniprot ID</b>	P15391
<b>Description</b>	PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Sterility</b>	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 $\mu$ m) prior to use.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only






Figure 1. Flow cytometry analysis with 1µl/test PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb on Raji cells (Blue histogram) or 293T cells (Red histogram).

