Cat. No. DMC100280P



PRODUCT INFORMATION

Clone ID **DMC280** CD21 **Target**

Synonyms C3DR; CD21; CR; CVID7; SLEB9

Host Species

PE-conjugated Anti-CD21 antibody(DMC280); Description

IgG1 Chimeric mAb

Delivery Under Development

Uniprot ID P20023

Rabbit/Human Fc chimeric IgG1 IgG type

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Reconstitution

Background

DIMA Disclaimer

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

> This gene encodes a membrane protein; which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been

found for this gene.

Usage Research use only

Conjugate PE-conjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are

actively scrutinizing all patent application to

ensure no IP infringement.







