

PRODUCT INFORMATION

Common Name	Iomab-B
Conjugate	PE-conjugated
Synonyms	B220;CD45;CD45R;GP180;L-CA;LCA;LY5;T200
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Mus musculus
IgG type	Mouse IgG1 - Kappa
Reactivity	Human
Target	CD45
Uniprot ID	P08575
Description	PE-conjugated Anti-CD45 (apamistamab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

