

PRODUCT INFORMATION

Clone ID DM118 CD7 **Target**

Synonyms CD7;GP40;TP41;LEU-9;Tp40

Host Species Rabbit

PE-conjugated Anti-CD7 antibody(DM118); Rabbit Description

mAb

Delivery Under Development

Uniprot ID P09564 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended Flow Cyt 1:100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Liquid PBS with 0.05% Proclin300, 1% BSA Reconstitution

Storage & Shipping Store at 2°C-8°C for 6 months

> This gene encodes a transmembrane protein which is a member of the immunoglobulin

superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell:B-cell **Background**

interaction during early lymphoid development.

Usage Research use only Conjugate PE-conjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under

Email: info@dimabio.com Website: www.dimabio.com

patent application. Any protein sequencing or **DIMA Disclaimer** reverse engineering attempt is prohibited. We are

actively scrutinizing all patent application to

ensure no IP infringement.

