

PRODUCT INFORMATION

Common Name	SGN-70, h1F6
Conjugate	PE-conjugated
Synonyms	CD70;CD27L;CD27LG;TNFSF7
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	CD70
Uniprot ID	P32970
Description	PE-conjugated Anti-CD70 (vorsetuzumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

