

PRODUCT INFORMATION

Clone ID 1E9 CD72 **Target** Lyb-2 **Synonyms Host Species** Rabbit

PE-conjugated Anti-CD72 antibody(1E9), IgG1 Description

Chimeric mAb

Delivery Under Development

Uniprot ID P21854

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended Flow Cyt 1:100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Liquid PBS with 0.05% Proclin300, 1% BSA Reconstitution

Storage & Shipping Store at 2°C-8°C for 6 months

Plays a role in B-cell proliferation and **Background**

differentiation.[UniProtKB/Swiss-Prot Function]

Usage Research use only

Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under

patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. **DIMA Disclaimer**

Email: info@dimabio.com Website: www.dimabio.com



