

**PRODUCT INFORMATION**

<b>Common Name</b>	BMS 986179
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	NT5E, 5'-Nucleotidase, 5'-NT, NT5, NTE
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Homo sapiens
<b>IgG type</b>	Human IgG2(C223S,Q274K,F296Y,F300Y,V309L,G327A,A330S,P331S,T339A,S378A,M397V) - Kappa
<b>Reactivity</b>	Human
<b>Target</b>	CD73
<b>Uniprot ID</b>	P21589
<b>Description</b>	PE-conjugated Anti-CD73(BMS 986179) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

