

PRODUCT INFORMATION

Common Name	CDH17-ADC-653-h43
Conjugate	PE-conjugated
Synonyms	Cadherin-17
Applications	Flow Cyt
Recommended Dilutions	
Formulation & Reconstitution	Liquid□PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	CDH17
Uniprot ID	Q12864
Description	PE-conjugated Anti-CDH17(TORL-3-600 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1  $\mu$ l/test PE-conjugated Anti-CDH17(TORL-3-600 biosimilar) on SNU-5 cell line (Blue histogram) or 293T cell line (Red histogram).

