

PRODUCT INFORMATION

Common Name	BIIB 059, BIIB-059, BIIB059
Conjugate	PE-conjugated
Synonyms	BDCA-2;CD303
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1(E356D,M358L) - kappa
Reactivity	Human
Target	CLEC4C
Uniprot ID	Q8WTT0
Description	PE-conjugated Anti-CLEC4C (litifilimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

