

PRODUCT INFORMATION

Common Name	5A11
Conjugate	PE-conjugated
Synonyms	CD355
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1(K214R,L234A,L235A,E356D,M358L) - Kappa
Reactivity	Human
Target	CRTAM
Uniprot ID	O95727
Description	PE-conjugated Anti-CRTAM(Oxford Biotherapeutics 5A11) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

