Cat. No. DMC100466P



PRODUCT INFORMATION

Clone ID DMC466 Target CSPG4

CSPG4A; HMW-MAA; MCSP; MCSPG; MEL-CSPG; **Synonyms**

MSK16; NG2

Host Species

PE-conjugated Anti-CSPG4 antibody(DMC466); IgG1 Chimeric mAb **Description**

Delivery Under Development

Uniprot ID O6UVK1

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

Dilutions

Flow Cyt 1:100

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Reconstitution

Background

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

> A human melanoma-associated chondroitin sulfate proteoglycan plays a role in stabilizing cell-substratum interactions during early events of melanoma cell spreading on endothelial basement membranes. CSPG4 represents an integral membrane chondroitin sulfate

proteoglycan expressed by human malignant melanoma cells. [provided by RefSeq; Jul 2008]

Usage Research use only Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under

patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to **DIMA Disclaimer**

Email: info@dimabio.com Website: www.dimabio.com

ensure no IP infringement.

