

PRODUCT INFORMATION

Common Name	Fab C225,IMC-225
Conjugate	PE-conjugated
Synonyms	EGFR;ERBB;ERBB1;HER1;PIG61;mENA
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 0.5 μ l/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Chimeric
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	EGFR
Uniprot ID	P00533
Description	PE-conjugated Anti-EGFR (Cetuximab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 0.25 μ l/test PE-conjugated Anti-EGFR (Cetuximab biosimilar) mAb on SNU-5 cells (Blue histogram) or Jurkat cells (Red histogram).

