

PRODUCT INFORMATION

Common Name	MFGR1877S
Conjugate	PE-conjugated
Synonyms	FGFR-3; CD333; JTK4
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	FGFR3
Uniprot ID	P22607
Description	PE-conjugated Anti-FGFR3(vofatamab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

