

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | OTSA101, ONC-005 |
| Conjugate | PE-conjugated |
| Synonyms | CD350, FZ-10, Fz10, FzE7, hFz10 |
| Applications | Flow Cyt |
| Endotoxin | Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization. |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Chimeric |
| IgG type | Human IgG1(E356D,M358L) - kappa |
| Reactivity | Human |
| Target | FZD10 |
| Uniprot ID | Q9ULW2 |
| Description | PE-conjugated Anti-FZD10(tabituximab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Sterility | Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

