

PRODUCT INFORMATION

Common Name	CAM-3001
Conjugate	PE-conjugated
Synonyms	GM-CSF-R-alpha; GMCSFR-alpha; GMR-alpha; CDw116; CD116; CSF2RA; CSF2R; CSF2RY
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG4(S228P) - Kappa
Reactivity	Human
Target	GM-CSFR
Uniprot ID	P15509
Description	PE-conjugated Anti-GM-CSFR(mavrilimumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

