Cat. No. BME100196P



PRODUCT INFORMATION

Common Name LBL-034

Conjugate PE-conjugated GPRC5D; CD3 **Synonyms Applications** Flow Cyt

Recommended Flow Cyt 1:100 **Dilutions**

Formulation & Liquid PBS with 0.05% Proclin300, 1% BSA

Reconstitution **Host Species** Humanized

lgG1 **IgG** type Reactivity Human

Target GPRC5D and CD3 Q9NZD1 and P07766 **Uniprot ID**

PE-conjugated Anti-GPRC5D and CD3(LBL-034 **Description**

biosimilar) BsAb

Delivery Under Development

Storage & Shipping Store at 2°C-8°C for 6 months

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

or animals.

Usage Research use only

