

PRODUCT INFORMATION

Common Name	LBL-034
Conjugate	PE-conjugated
Synonyms	GPRC5D; CD3
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	GPRC5D and CD3
Uniprot ID	Q9NZD1 and P07766
Description	PE-conjugated Anti-GPRC5D and CD3(LBL-034 biosimilar) BsAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

