

PRODUCT INFORMATION

Common Name	Innovative Cellular Therapeutics
Conjugate	PE-conjugated
Synonyms	DIAR6, GC-C, GUC2C, MECIL, MUCIL, STAR
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	GUCY2C
Uniprot ID	P25092
Description	PE-conjugated Anti-GUCY2C (ICTCAR-CRC) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

