

PRODUCT INFORMATION

Common Name	AIN457
Conjugate	PE-conjugated
Synonyms	CTLA-8;CTLA8;IL-17;IL-17A;IL17
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL17A
Uniprot ID	Q16552
Description	PE-conjugated Anti-IL17A(secukinumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

