

PRODUCT INFORMATION

Common Name	CAN-04, CAN04
Conjugate	PE-conjugated
Synonyms	IL1AP; C3orf13; IL1R3; IL-1R-3; IL-1R3
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL1RAP
Uniprot ID	Q9NPH3
Description	PE-conjugated Anti-IL1RAP(nadunolimab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 100 μ l/test (1:400) PE-conjugated Anti-IL1RAP(nadunolimab biosimilar) mAb (BME100241P) on HEK293 cells transfected with Human IL1RAP and IL1RA protein (Blue histogram) or HEK293 transfected with Human IL1RA protein (Red histogram).

