

PRODUCT INFORMATION

Common Name	BI-655130, BI655130, BI 655130, spesolimab-sbzo
Conjugate	PE-conjugated
Synonyms	ILRL2; IL-36R; IL1RRP2; IL-36R
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	Human IgG1(K214R, L234A, L235A) - Kappa
Reactivity	Human
Target	IL1RL2
Uniprot ID	Q9HB29
Description	PE-conjugated Anti-IL1RL2(spesolimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

