

PRODUCT INFORMATION

Common Name	BI-655130,BI655130,BI 655130,spesolimab-sbzo
Conjugate	PE-conjugated
Synonyms	ILRL2; IL-36R; IL1RRP2; IL-36R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1(K97R,L117A,L118A)
Reactivity	Human
Target	IL1RL2
Uniprot ID	Q9HB29
Description	PE-conjugated Anti-IL1RL2(spesolimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

