

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Common Name                  | 15D2, NN-8226, NNC-0109-0012  |
| Conjugate                    | PE-conjugated   |
| Synonyms                     | IL-20;Cytokine Zcyto10  |
| Applications                 | Flow Cyt  |
| Recommended Dilutions        | Flow Cyt 1:100  |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA  |
| Host Species                 | Homo sapiens  |
| IgG type                     | IgG4  |
| Reactivity                   | Human   |
| Target                       | IL20  |
| Uniprot ID                   | Q9NYY1  |
| Description                  | PE-conjugated Anti-IL20(fletikumab biosimilar) mAb  |
| Delivery                     | Under Development   |
| Storage & Shipping           | Store at 2°C-8°C for 6 months   |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage                        | Research use only   |

