

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | 15D2, NN-8226, NNC-0109-0012 |
| Conjugate | PE-conjugated |
| Synonyms | IL-20;Cytokine Zcyto10 |
| Applications | Flow Cyt |
| Endotoxin | Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization. |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | Human IgG4 - Kappa |
| Reactivity | Human |
| Target | IL20 |
| Uniprot ID | Q9NYY1 |
| Description | PE-conjugated Anti-IL20(fletikumab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

